

Model





CABINET **MOUNTED** SYSTEMS

Flexible options Specialty solutions Cabinet mounted hand controls



CABINET MOUNTED SYSTEMS

With an industry-leading assortment of positioning options, BDS specializes in cabinet mounted handpiece controls. Whether the requirement is rear delivery, side delivery, arm mounted, or shelf mounted, BDS has every option needed to provide specialty solutions for every practice and preference.

Model SC-4250

- Doctor's side = Three-handpiece automatic with air/water controls
- Assistant's side = Three-way syringe, autoclavable S/E and HVE valves
- System Purge feature flushes entire water system and replaces with clean water
- Mounts to underside of counter top with included hardware
- Easy left-right positioning

Model SC-4200

- Three-handpiece automatic with air/water controls
- System Purge feature flushes entire water system and replaces with clean water
- Mounts to underside of counter top with included hardware
- Side mount, recessed controls

Model **A-5150**

- Three-way syringe, autoclavable S/E and HVE valves
- Mounts to underside of counter top with included hardware







Model **S0-3304, S0-3354** (shown)

- Three-handpiece automatic with air/water controls
- Autoclavable S/E and HVE valves (SO-3354)
- System Purge flushes entire system & replaces with clean water
- System extends 12" from recessed position & pivots 30° left-right
- Mounts to underside of counter top with included hardware
- Side mount, recessed controls
- SO-3304 = Without HVE and S/E valves
- SO-3354 = With HVE and S/E valves

Model SD-3302

- Three-handpiece automatic with air/water controls
- System Purge feature flushes entire water system and replaces with clean water
- Heavy-duty, pneumatic air-brake flex-arm (9" or 24" rigid and 24" flex section with 25" vertical travel)
- Side mount, recessed controls
- Mounts to side or rear of cabinet with included hardware

Model

A-230, A-250 (shown)

- Two-handpiece with air/water controls (one wet and one dry)
- Automatic handpiece selection (A-250)
- Designed for under-shelf mounting
- A-230 = Without automatic handpiece selection
- A-250 = Automatic handpiece selection







FEATURES OF CABINET MOUNTED SYSTEMS	SC-4250	SC-4200	A-5150	A-230	A-250	SO-3304	SO-3354	SD-3302
2 handpiece automatic (1 wet, 1 dry)Non-retracting 2 2								
2 handpiece manual (1 wet, 1 dry)								
3 handpiece automatic								
3-way air/water syringe								
Master on/off toggle with indicator								
Pressure adjustment for each handpiece								
Easy to read pressure gauge								
Non-retracting water coolant, adjustment & purge								
Non-retracting water coolant & adjustment								
Air coolant adjustment								
Asepsis handpiece tubing								
Vacuum canister with removable screen (131-052)								
Vacuum canister with removable screen (131-053)								
1 Autoclavable HVE and saliva ejector valve w/ tubing								¢
Asepsis flex-arm system with pneumatic brake								
2″ dia. wall mount bracket								
Slide-out mount with 15° pivot left and right								
Standard disc foot control								
Master valve assembly with filter and regulators								
7' umbilical tubing								

Please contact BDS for technical cut-out templates.



800.237.2303

www.bdsystems.us

800.237.2303 info@beaverstatedental.com