





EPICOPERATORY CHAIR



EPICOPERATORY CHAIRS

Power. Strength. Durability. The Epic Chair is the high-water mark of toughness. Our chair leads the way in dependability and does it with the style and comfort that will keep both staff and patients happy. Guaranteed to lift 600 pounds



over and over again, the Epic Chair truly raises the bar. With a 7-year warranty, and a 10-year warranty on the electromechanical motors, BDS provides a chair that will be there for the long haul.

Epic Chair options

- Ultraleather upholstery upgrade
- Additional touchpad
- Foot control
- Replacement toe covers

Chair base floor pad CH-0893




Foot switch with additional preset **400-045**



Second touch-pad armrest control **400-165**





FEATURES OF EPIC OPERATORY CHAIR	
Heavy duty cast iron construction	_
This back with lumbar support	
Armrest mounted touch pad	
4 manual chair movements, 3 programmable positions & home position	
Swing down armrest	
Synchronized back and seat recline	
Heavy vinyl scuff guard for toe board	
Electromechanical chair operation (110 VAC)	
600 lb. lifting capacity	
Double articulating headrest	
30° Rotation left and right of center	
Standard Naugasoft Upholstery Options (see BDS catalog page 128)	
OPTIONS OF EPIC OPERATORY CHAIR	EPIC CHAIR
Epic chair foot control	
Additional touch pad	
Ultraleather™ (in lieu of standard upholstery)	

Please contact BDS for technical cut-out templates.

