



AMERICAN MADE

# BDS<sup>®</sup>



## HANDPIECE CONTROL SYSTEMS

Simple & compact controls

Cost effective solutions

Versatile for all practices & labs



Model  
S-5325



Model  
S-330

# HANDPIECE CONTROL SYSTEMS

When an operatory is down, but a completely new package is not needed, BDS has retrofitting options to get that operatory back up and running quickly, efficiently, and economically. From single handpiece

## Model S-330

- Three-handpiece automatic with air/water controls
- Rugged aluminum cover
- Side mount, recessed controls
- Non-retracting water coolant and purge



*Unit head dimensions (w x d) 14" x 13"*

## Model S-5325

- Two-handpiece automatic (1 wet, 1 dry)
- Easy to read pressure gauge
- Vacuum canister with removable screen (131-053)
- Telescoping arm with breakaway feature, 29" to 40" horizontal reach
- Non-retracting water coolant



*Unit head dimensions (w x d) (11" x 6½")*

## Model A-300

- Three-handpiece automatic with air/water controls
- Pressure adjustment for each handpiece
- Non-retracting water coolant and purge



operation to a three-handpiece control with an air/water syringe, BDS's wide range of simple and compact controls will fill the need of any dental practice or lab.

### Model

## A-220, A-222, A-200, A-202

- Two-handpiece automatic
- Pressure adjustment for each handpiece
- Non-retracting water coolant for 2 handpieces A-222
- Non-retracting water coolant for 1 handpiece (1 wet, 1 dry) A-220
- A-220 = 2 Handpiece automatic
- A-222 = 2 Handpiece automatic both wet (not shown)
- A-200 = 2 Handpiece manual
- A-202 = 2 Handpiece manual both wet (not shown)



Dimensions (w x h x d)  
6 1/8" x 5 3/4" x 3 1/2"



Unit head dimensions (w x h x d)  
7" x 5 3/4" x 3 1/2"

### Model A-400

- Three-handpiece automatic
- Pressure adjustment for each handpiece



Unit head dimensions (w x h x d)  
6 1/8" x 5 3/4" x 3 1/2"


### Model A-100

- Single handpiece control with coolant
- Handpiece pressure adjustment
- Non-retracting water coolant



### Model A-101

- Single handpiece control with out coolant
- Handpiece pressure adjustment
- Easy to read pressure gauge



### Model A-050

- Single handpiece control with out coolant
- Handpiece pressure adjustment
- Easy to read pressure gauge

800.237.2303

info@beaverstatedental.com



[www.bdsystems.us](http://www.bdsystems.us)  
115 S. Elliott Road, Newberg, OR 97132  
800.237.2303

800.237.2303  
[info@beaverstatedental.com](mailto:info@beaverstatedental.com)