



# BDS



## OVER THE PATIENT SYSTEMS

Clean and contemporary  
Wide range of configurations  
New operator to replacement equipment



Model  
S-3600

# OVER-THE-PATIENT SYSTEMS

With a wide variety of arm-mounted systems, BDS has the Over-the-Patient solutions you need. From a full new operatory package with traditional three handpiece

## Model

### S-3600, S-3606

- Three-handpiece automatic with air/water controls
- Plastic cover with smooth edges
- Ergonomic, easy-access controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water
- 07-145 Unit mount touch pad available on S-3600, S-3606, S-3630 & S-3636 only.

■ S-3600 = Over-post mount (shown) ■ S-3606 = Top-post mount



Unit head dimensions (w x d) 15.5" x 15.5"

## Model

### S-3630, S-3636

- Three-handpiece automatic with air/water controls
- Plastic cover with smooth edges
- Ergonomic, easy-access controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water
- 07-145 Unit mount touch pad available on S-3600, S-3606, S-3630 & S-3636 only.
- Unit mount swivel tray 9-3/4" x 13-1/2"

■ S-3630 = Over-post mount (shown) ■ S-3636 = Top-post mount




Unit head dimensions (w x d) 15.5" x 15.5"

## Model


### S-3420

- Three-handpiece automatic with air/water controls
- Manual counter-balanced flex-arm (27.5" rigid and 26.5" flex sections with 27" vertical travel)
- Suitable for doctor or hygienist use

■ S-3420 = Over-post mount



Unit head dimensions (w x d) 10.5" x 8"



## Unit Mount Touch Pad 107-145

- available on S-3600, S-3606, S-3630 & S-3636

system, to retro-fitting an existing chair with a European style system to simply hygiene controls, BDS gets operatories running quickly, efficiently, and reliably.

## Model

### S-6300, S-6307

- Three-handpiece automatic with air/water controls
- Ergonomic and functional Euro styling
- Side mount, recessed controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water

■ S-6300 = Over-post mount ■ S-6307 = Top-post mount (shown)



Unit head dimensions (w x d) 13.25" x 11.5"

## Model

### S-3300, S-3306

- Three-handpiece automatic with air/water controls
- Rugged aluminum cover
- Side mount, recessed controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water

■ S-3300 = Over-post mount (shown) ■ S-3306 = Top-post mount




Unit head dimensions (w x d) 14" x 14"

## Model

### S-2500

- Two-handpiece automatic with air/water controls
- Manual counter-balanced flex-arm (27.5" rigid and 26.5" flex sections with 27" vertical travel)
- Fixed instrument tray
- Designed for hygienist use

■ S-2500 = Over-post mount



Unit head dimensions (w x d) 8" x 6"

800.237.2303

info@beaverstatedental.com

FEATURES OF OVER-THE-PATIENT SYSTEMS	S-3600	S-3606	S-3630	S-3636	S-3300	S-3306	S-6300	S-6307	S-3420	S-2500
2 handpiece automatic control										■
3 handpiece automatic control	■	■	■	■	■	■	■	■	■	
3-way air/water syringe	■	■	■	■	■	■	■	■	■	■
Master on/off toggle with indicator	■	■	■	■	■	■	■	■		
Master on/off toggle									■	■
Pressure adjustment for each handpiece	■	■	■	■	■	■	■	■	■	■
Easy to read pressure gauge	■	■	■	■	■	■	■	■	■	■
Non-retracting water coolant with adjustment									■	■
Non-retracting water coolant adjustment and purge	■	■	■	■	■	■	■	■		
Air coolant adjustment	■	■	■	■	■	■	■	■		
Asepsis handpiece tubing	■	■	■	■	■	■	■	■	■	■
Stainless Steel instrument tray, with non-slip pad					■	■	■	■		
Asepsis flex-arm system with pneumatic brake	■	■	■	■	■	■	■	■		
Manual flex-arm system with friction brake									■	■
Unit mount post, 28"		■		■		■		■		
Standard disc foot control	■	■	■	■	■	■	■	■	■	■
Utility center assembly with master valves & regulators	■	■	■	■	■	■	■	■	■	■
7' umbilical tubing	■	■	■	■	■	■	■	■	■	■

Please contact BDS for technical cut-out templates.



www.bdsystems.us  
 115 S. Elliott Road, Newberg, OR 97132  
 800.237.2303

800.237.2303

info@beaverstatedental.com