





REAR CABINET MOUNTED SYSTEMS

Compact and versatile
Cost effective
Configure to positioning preferences



REAR CABINET MOUNTED SYSTEMS

If you are looking for the versatility of a rear cabinet mounted system, BDS offers a wide array of options for you. Configure it how you like and put it to work. BDS systems will keep your operatory working reliably, efficiently, and effectively.

Model

RD-4600, RD-4650 (shown)

- Under counter mount duo system (RD-4650)
- Under counter mount doctor's system (RD-4600)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-4650)



Model

RD-4250, **RD-4251** (shown)

- Wall/cabinet mount duo system (RD-4250)
- Wall/cabinet mount Assistant's system (RD-4251)
- Three handpiece automatic with air/water controls (RD-4250)
- Autoclavable SE and HVE
- Starboard work surface
- Industrial strength pneumatic support arm
- Left/right capability (RD-4250)



Model

RD-3103, **RD-3153** (shown)

- Under counter mount duo system (RD-3153)
- Under counter mount doctor's system (RD-3103)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-3153)



Model

RD-3100, RD-3150 (shown)

- Under counter mount duo system with left/right swing
- Three handpiece automatic with air/water controls
- Autoclavable SE and HVE
- Left/right capability
- Standard disc foot control
- Air coolant adjustment
- RD-3100 = No vacuum RD-3150 = With vacuum



Model

RD-3101, RD-3151 (shown)

- Under counter mount duo system (RD-3151)
- Under counter mount doctor's system (RD-3101)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-3151)



800.237.2303

info@beaverstatedental.com

FEATURES OF REAR CABINET MOUNTED SYSTEMS	RD-3150	RD-3100	RD-3101	RD-3151	RD-3103	RD-3153	RD-4250	RD-4251	RD-4600	RD-4650
3 handpiece automatic										
Non-retracting water coolant adjustment and purge										
Pressure adjustment for each handpiece										
Dual 3-way air/water syringes										
3-way air/water syringe										
Master on/off toggle with indicator										
Asepsis handpiece tubing										
Autoclavable, HVE saliva ejector valve & tubing										
Wall mounting plate										
Easy to read pressure gauge										
Standard disc foot control										
Vacuum canister with removable screen (131-052)										
Vacuum canister with removable screen (131-053)										
Master valve assembly with filter and regulators										
7' umbilical tubing										

Please contact BDS for technical cut-out templates.

